

FIGURE 7. LFA-1 biosynthesis in patient or health adult lymphocytes. PHA-stimulated lines (A, B) or EBV-transformed lines (C-E) were labeled with [55 S]methionine for 2 h, and either harvested immediately or chased for 22 h with unlabeled methionine as indicated. Lysates with Triton X-100 and sodium deoxycholate were immunoprecipitated with purified antibodies coupled to Sepharose CL-4B: TS/22 anti-LFA-1 α L, TS1/18 anti- β , or activated, quenched Sepharose, which served as a negative control. Precipitates were subjected to SDS-6% PAGE and fluorography. Only the upper part of the gel is shown.