In This Issue

Protrusion provides a prognosis

ow a cancer cell migrates across a 2D surface says little about its ability to move through a 3D matrix, Meyer et al. reveal, but measuring membrane protrusion can be much more informative.

A metastasizing tumor cell moves through 3D tissues, yet, largely because they're easier to perform, most studies of cell migration are carried out on 2D surfaces. Meyer et al. wanted to investigate how well 2D assays replicate 3D motility and to determine what parameters best predict a cell's migration through 3D environments.

The researchers tracked the migration of several breast cancer cell lines on a variety of 2D and 3D matrices in response to different growth factors associated with metastasis. The speed and directional persistence of the cells' movements on 2D surfaces bore little resemblance to their behavior in 3D matrices. And, although a cell's

migratory response to different growth factors could be partially predicted by measuring which growth factor receptors it expressed, differences in receptor expression and activation levels couldn't account for all the variation seen in different cell lines.

Meyer et al. then examined whether any individual steps of 2D migration correlated with 3D motility and found that the extent to which cells formed membrane protrusions in response to different growth factors was a reliable predictor of their 3D migratory behavior. Strengthening the connection, cytoskeletal inhibitors that preferentially block 3D migration also impeded membrane protrusion in 2D. The authors now want to develop high-throughput methods for measuring membrane protrusion in order to facilitate screens for drugs and siRNAs that affect 3D cell migration.

Meyer, A.S., et al. 2012. *J. Cell Biol.* http://dx.doi.org/10.1083/jcb.201201003.

The ribosome runs a systems test



Within the ribosomal large subunit, Rpl10 (blue) embraces the P-site tRNA (purple). If the P-site is intact, Efl1 (red) induces the release of Tif6 (green).

B ussiere et al. describe how new ribosomes check that they're ready to start synthesizing proteins.

The large and small subunits of

eukaryotic ribosomes are preassembled in the nucleus, but, after their export into the cytoplasm, they undergo a series of maturation steps before they start translating mRNAs. In one of the final maturation events, the large subunit sheds an inhibitory protein called Tif6 so

that it can associate with the small sub-

unit and initiate protein synthesis. Bussiere et al. found that Tif6 wasn't released from yeast large subunits lacking the ribosomal protein Rpl10.

Rpl10 is located next to the P-site—the catalytic heart of ribosomes where amino acid-carrying tRNAs are coupled to the grow-

ing polypeptide chain. Mutating a region of Rpl10 that contacts the P-site tRNA also blocked Tif6 release, thereby inhibiting yeast cell growth. Maturation and growth were restored by mutations in Tif6 that weaken the inhibitory protein's association with the large subunit. Rpl10 mutants were also rescued by mutants of Efl1, a GTPase that promotes Tif6 release. The mutations in Efl1 appear to facilitate a conformational change analogous to changes observed in the translocation factor eEF2 during translation.

The results suggest that Efl1 normally removes Tif6 from large subunits only if the P-site is correctly assembled, ensuring that new ribosomes don't become active if their catalytic sites are nonfunctional. Senior author Arlen Johnson now wants to examine the role of Sdo1, another protein required for Tif6 release that is shaped somewhat like a tRNA and that may therefore help Efl1 to determine whether or not the P-site is able to support protein synthesis.

Bussiere, C., et al. 2012. *J. Cell Biol.* http://dx.doi.org/10.1083/jcb.201112131.

IFT-A goes both ways





The primary cilium looks fairly normal in a weak ift144 mutant (left) but is greatly truncated when the protein is absent (right).

he intraflagellar transport A (IFT-A) complex regulates Sonic hedgehog (Shh) signaling by trafficking membrane proteins into primary cilia, Liem et al. reveal.

Vertebrate cells require a primary cilium to activate Shh signaling because many of the pathway's key components accumulate inside these mi-

crotubule-based organelles. Primary cilia themselves rely on two multi-subunit transport complexes, IFT-A and IFT-B. IFT-B moves proteins to the tips of cilia, and, in the absence of the complex, cells are unable to form cilia or activate the Shh pathway. IFT-A's main role, on the other hand, has been thought to be the recycling of proteins out of the primary cilium. Mice lacking certain IFT-A subunits have swollen primary cilia and show ectopic Shh activity.

Liem et al. identified two mouse lines with different mutations in

the IFT-A subunit *ift144*. One line, named *twinkletoes*, expressed a weakly active version of IFT144 and was similar to other IFT-A mutant mice in having slightly larger primary cilia and enhanced Shh signaling. The *diamondhead* line, however, produced no IFT144 protein, had very short cilia, and showed a decrease in Shh activity. This suggests that, in addition to its recycling function, IFT-A transports proteins into cilia and, like IFT-B, is required for cilia assembly. Indeed, several membrane-associated components of the Shh pathway were missing from the short cilia of *diamondhead* mice. Soluble pathway components were still present in cilia, however, perhaps indicating a division of labor between IFT-A and IFT-B.

Shh activity may be increased in *twinkletoes* mice because less adenylyl cyclase is transported into their primary cilia. The researchers propose that this membrane-associated enzyme normally generates cyclic AMP in order to activate protein kinase A, a potent inhibitor of the Shh pathway.

Liem, K.F., Jr., et al. 2012. *J. Cell Biol.* http://dx.doi.org/10.1083/jcb.201110049.